

*Instructions for SOTA Bio Tuner  
Model BT6pro*



**CONTENTS**

WELCOME TO SOTA ..... 2

BEFORE USING ..... 2

LEARN MORE ..... 2

COMPLETE UNIT INCLUDES..... 2

DO NOT USE..... 3

CAUTIONS ..... 3

SUMMARY OF LIGHTS ..... 4

HOW TO USE THE BIO TUNER..... 4

    Headset Option..... 4

    Conductive Pads Option..... 5

    Ear Clips Option..... 5

    Turning the Unit On ..... 5

    How it Should Feel..... 7

    Typical Use..... 7

OUTPUT OF SIX MODES..... 8

ELECTRICAL SPECIFICATIONS ..... 9

    Power Requirements..... 9

    Optional Power Input..... 9

    Output Specifications..... 9

REPLACEMENT ACCESSORIES ..... 9

WARRANTY..... 9

REPAIRS..... 9

• The SOTA units are consumer products designed to be used as part of a Wellness Lifestyle • The SOTA units are not medical devices nor are they intended to diagnose, treat or cure any medical or health condition • The use, safety or effectiveness of the SOTA units has not been approved by any government agency • Please consult a health practitioner for any medical or health condition. • Results using the SOTA Products will vary for each person. • Please note that all purchases of SOTA Products are deemed to be made in Penticton, BC, Canada •





## WELCOME TO SOTA

Thank you for allowing us to be a part of your wellness team.

The Bio Tuner is a consumer product designed to be used as part of a Wellness Lifestyle. A commitment to healthy lifestyle choices is an important factor in the journey to Wellness.

The SOTA Bio Tuner supports the mind-body connection. It outputs a broad range of harmonic frequencies helping to create an inner sense of calm, balance and tranquility.

The Bio Tuner is designed to be used independently, but can also be used together with any of the SOTA Products.

With proper care, your Bio Tuner should provide you many years of trouble-free use.

Please understand that results will vary.

## BEFORE USING

Please read this manual thoroughly before using your Bio Tuner.

Please note information in this manual is subject to change.

For the most up-to-date manual, please visit our website [www.sota.com/manuals](http://www.sota.com/manuals).

## LEARN MORE

To learn more about the Bio Tuner we invite you to visit the SOTA website:

[www.sota.com](http://www.sota.com)

The website offers How to Videos, Frequently Asked Questions, *The SOTA Products User Guide*, SOTA News ... and more.

## COMPLETE UNIT INCLUDES

- One (1) Bio Tuner Unit
- One (1) 9-Volt Alkaline Battery
- One (1) Product Manual





## HEADSET OPTION

- One (1) Headset
- Six (6) Cotton Sleeves

## CONDUCTIVE PADS OPTION

- One (1) Cord
- Ten (10) pair Disposable Conductive Pads

## EAR CLIPS OPTION

- One (1) Cord
- One (1) Pair of Ear Clips

## DO NOT USE

1. If you are prone to or experience tinnitus.
2. In the following situations as safety has not been explored:
  - a. In the case of pregnancy.
  - b. With an implanted device such as a pacemaker, etc.
3. Discontinue use if:
  - a. You experience headaches or pain when using the unit.
  - b. You experience a ringing in your ears when using the unit.
  - c. You experience an increase in neuropathy-related nerve pain when using the unit.
4. As with all electrical products do not use near water, when driving a car or when operating heavy equipment.
5. Do not:
  - a. Suddenly increase the intensity to the maximum level. The intensity should be increased gradually to a comfortable level to ensure safety.

## CAUTIONS

1. If you are on psychiatric medication or have a neurological condition, please consult your health care provider prior to using the Bio Tuner.
2. In very rare instances, individuals have reported the onset of tinnitus when using the unit. If you experience tinnitus or a ringing in the ears, please discontinue use.



3. In very rare instances, individuals have reported experiencing nerve pain when using the unit. If you experience nerve pain, please discontinue use.
4. While this technology is generally considered safe, there exists potential for rare individual reactions that cannot reasonably be foreseen. Therefore, your use of the SOTA Products constitutes your agreement that you are responsible for your decision to use the technology.

## SUMMARY OF LIGHTS


**GREEN:** Indicates the Bio Tuner is in operation and the current is flowing through the Headset. Also indicates Modes 1 through 3.

**RED:** Indicates the battery is low and will soon need replacing. Also indicates Modes 4 through 6.


## HOW TO USE THE BIO TUNER

1. Place the 9-Volt Alkaline battery in the unit.
  - a. On the back of the unit, slide off the battery cover.
  - b. Clip the battery connections in the battery compartment to the 9-Volt Alkaline battery terminals.
  - c. Slide the 9-Volt Alkaline battery into the case and push down. Slide the cover back onto the unit.

## HEADSET OPTION



2. Clean the area behind each ear to eliminate skin oils. This will allow for better conductivity.
  - a. Wet the cotton sleeves with water and place them over the Headset electrodes.
  - b. Insert the Headset Cord into the jack on the unit.
  - c. Position the Headset electrodes in the hollow behind the ears as in image 1.









## CONDUCTIVE PADS OPTION

2. Clean the area behind each ear to eliminate skin oils. This will allow for better conductivity.
  - a. Insert the Cord into the jack on the unit.
  - b. Attach one set of Conductive Pads to each connection on the Cord.
  - c. Attach a Conductive Pad behind each ear lobe as in image 2.



## EAR CLIPS OPTION

2. Clean the ear lobes to eliminate skin oils. This will allow for better conductivity.
  - a. Insert the pins on the cord into the Ear Clips.
  - b. Wet the conductive rubber on the Ear Clips with water.
  - c. Place one Ear Clip on each ear lobe, positioning the Ear Clip with the conductive rubber at the back of the ear lobe and close to the jaw as in image 3.



## TURNING THE UNIT ON

3. Turn the unit ON by turning the ON/OFF switch clockwise.

The GREEN light on the end panel will begin to flash alternately once every second, indicating the unit is in Mode 1. If the lights do not come on, the battery may not be properly connected.

4. Gradually turn up the intensity to a comfortable level.  
The intensity should be comfortable. See 'How It Should Feel' on page 7.
5. Place the unit in a pocket.
6. Change Modes by pushing the Mode Select button. The table on page 6 indicates how to identify each Mode. For more information on the Modes, see Output of Six Modes on page 8.



| MODE   | INDICATED BY                                    |
|--------|---|
| MODE 1 | GREEN light flashes once (1) per second.        |
| MODE 2 | GREEN light flashes twice (2) times per second. |
| MODE 3 | GREEN light flashes three (3) times per second. |
| MODE 4 | RED light flashes once (1) per second.          |
| MODE 5 | RED light flashes twice (2) times per second.   |
| MODE 6 | RED light flashes three (3) times per second.   |

7. After the 20-minute timed session, the unit will automatically cycle into Sleep Mode and both the GREEN and the RED lights will stop flashing and come on steady.
8. While in Sleep Mode the battery is still being used.
  - a. To turn the unit OFF, rotate the ON/OFF switch counter clockwise until it clicks. Do not force the ON/OFF switch past its natural stopping point.
  - b. To continue with another session, turn the unit ON again and choose your desired Mode. Before starting a second timed cycle, wet the cotton sleeves again.

**Notes:**

**Headset and Ear Clips:** The Headset and the Ear Clips may feel more comfortable with the use of either conductive gel or aloe vera gel rather than water. Over time, the Ear Clips may lose their conductivity and will need to be replaced.

**Changing Modes:** When the Mode is changed during a timed session, the timing starts again for a new 20 minute session.

**Care of Bio Tuner Cord and Ear Clips:** Occasionally, remove the conductive rubber from the Ear Clips to clean in a mild soap solution. Rinse thoroughly. If necessary, gently wipe the metal pins of the Bio Tuner Cord.

**Low Battery Warning:** The RED light comes on when the battery is low. This means the battery will soon need to be replaced.



## HOW IT SHOULD FEEL

1. The sensation is a pleasant buzzing or the current may not be felt at all. It is not necessary to turn the control up to an uncomfortable level. If it feels prickly or unpleasant, check the following:

### HEADSET OPTION

- a. Wet the cotton sleeves again.
- b. Apply conductive gel or aloe vera gel to the hollow behind the ears.

### CONDUCTIVE PADS OPTION

- a. Clean the area behind the ears again.
- b. Replace the Conductive Pads.

### EAR CLIPS OPTION

- a. Wet the conductive rubber of the Ear Clips again.
  - b. Apply conductive gel or aloe vera gel to the ear lobes.
3. With repeated use, the buzzing sensation may not be felt as much. The intensity may be turned up if desired.
  4. You may feel a sensation on one side more than the other.
  5. In rare instances, some individuals may not feel the current.

## TYPICAL USE

**Length of Session:** One cycle is timed for approximately 20-minutes. One session can be 20 to 40 minutes (1 to 2 cycles).

Information indicates that if used longer than 40 minutes per session, the benefits may be reduced. The unit can be used for up to three 40-minute sessions per day – space the sessions throughout the day.

**Length of Use:** The unit should be used daily, at least 2 sessions per day, for a minimum of 4 weeks or longer. After the benefits have been achieved, simply use the unit as needed.

**Intensity Level:** The intensity should always be set to a comfortable level. Be sure to reduce the intensity before changing modes. After the new mode is selected, readjust the intensity.

**With Extended Use:** It is recommended to vary the Mode from time to time.



**Sensitive Skin:** During the first weeks of use, the skin may become irritated when using the Ear Clip Option. Over time, irritation usually decreases or disappears. Apply a healing lotion after use if the ear lobes become red or irritated.

**Drink plenty of water:** It's important to drink enough water when using the unit. Drinking ozonated water is especially helpful.

**Stop Gradually:** Reduce the amount of time each day for the last week of use.

## OUTPUT OF SIX MODES

The Bio Tuner provides six different Modes with which to experiment and discover what feels or works best. Try different Modes to determine which feels most comfortable. With extended use, it is recommended to vary the Mode from time to time, as your body can become acclimatized.

One Mode is not considered more beneficial than another.

In all Modes, the base frequency of 1,000 Hz pulse-modulated by 111.11 Hz is present. Modes 4, 5, and 6 may have a more gentle feel.

See table of Modes below.

### BIO TUNER MODES

#### MODE 1-3 (550 uS pulse width)

Mode 1: 1,000 Hz modulated by 111.11 Hz,  
Output = Continuous

Mode 2: 1,000 Hz modulated by 111.11 Hz,  
Output = 1 Second ON, 1 Second OFF (0.5 Hz)

Mode 3: 1,000 Hz modulated by 111.11 Hz,  
Output = Amplitude Modulated by 7.83 Hz Schumann (Earth) Sine Wave

#### MODE 4-6 (220 uS PULSE WIDTH)

Mode 4: 1,000 Hz modulated by 111.11 Hz,  
Output = Continuous

Mode 5: 1,000 Hz modulated by 111.11 Hz,  
Output = 1 Second ON, 1 Second OFF (0.5 Hz)

Mode 6: 1,000 Hz modulated by 111.11 Hz,  
Output = Amplitude Modulated by 7.83 Hz Schumann (Earth) Sine Wave





## **ELECTRICAL SPECIFICATIONS**

### **POWER REQUIREMENTS**

Use a 9-Volt Alkaline battery for proper operation.

### **OPTIONAL POWER INPUT**

The Bio Tuner can also be operated using NiMH rechargeable batteries. When purchasing rechargeable batteries, it is suggested to take the unit with you to ensure the batteries will fit the enclosure.

### **OUTPUT SPECIFICATIONS**

Maximum <1.5mA Peak ~40 Volts Peak-Peak into typical 4k ohm load. Minimum <~10uA, 10 Volt Peak-Peak at lowest setting. Output impedance ~20k ohm.

1,000 Hz (+/- 1%) pulse-modulated by 111.11 Hz (+/- 1%). Waveforms are Quartz Crystal-Controlled for precise frequency output. There is also a broad range of harmonic frequencies delivered in each pulse as observed on a frequency domain spectrum analyzer.

Pulses are Bi-Phasic (AC) with no direct current (DC) or offset bias. Output is transformer coupled so that no DC can appear at electrodes.

## **REPLACEMENT ACCESSORIES**

To replace your accessories, contact your supplier.

## **WARRANTY**

The Bio Tuner, BT6pro is covered by a three (3) year limited warranty. Please keep proof of purchase. Warranty is void if the failure is due to abuse or negligence. SOTA reserves the right to make changes to the Bio Tuner without obligation.

## **REPAIRS**

Please contact your supplier for return instructions. Or, for information about your closest repair depot, contact [repairs@sota.com](mailto:repairs@sota.com) or

Toll-Free: 1.800.224.0242 (Canada & USA) Tel: 1.250.770.2023

Designed, Engineered and Made in Canada.